



ALPHAVAC

MULTIPURPOSE MECHANICAL ASPIRATION





CONTROLLED REMOVAL. MAXIMUM VERSATILITY. RESTORED FLOW.


The AlphaVac System is an efficient way to remove clot from large vessels.

THE ALPHAVAC SYSTEM **ALLOWS FOR**


- + Large vessel removal from neck to knee
- + Confirmation of clot engagement
- + Easy navigation and control
- + Enhanced torqueability
- + Removal of clots with limited blood loss



Superior vena cava



Inferior vena cava




Iliofemoral Vein

FEATURES



CANNULA	F22 ²⁰	F22 ¹⁸⁰	F18 ⁸⁵
PART NUMBER	H787253000	H787253010	H787253020
ACCESS	26F	26F	22F
CANNULA SIZE	22F	22F	18F
LENGTH	77cm	77cm	105cm
FUNNEL SIZE	14mm	14mm	11mm
ANGLE	20°	180°	85°

Experience the
AlphaVac System
Alpha-Vac.com



ACCESS SITES

AlphaVac Risk Information

Indications for Use:

The F18⁸⁵ Cannula is indicated for:

- the non-surgical removal of thrombi or emboli from venous vasculature.
- aspiration of contrast media and other fluids from venous vasculature.

The Cannula is intended for use in the venous system. The Handle is indicated as a vacuum source for the AlphaVac Multipurpose Mechanical Aspiration System.

The AlphaVac F22 Cannula is indicated for:


- the non-surgical removal of thrombi or emboli from the vasculature
- aspiration of contrast media and other fluids from the vasculature

The Cannula is intended for use in the venous system. The Handle is indicated as a vacuum source for the AlphaVac Multipurpose Mechanical Aspiration System

Contraindications: The device is contraindicated in the removal of chronic firmly adherent intravascular material (e.g., atherosclerotic plaque, chronic pulmonary embolism). The device is contraindicated for use in the right heart during active cardiopulmonary resuscitation. The device is contraindicated for blood storage and infusion back into the patient.

Refer to Directions for Use and/or User Manual provided with the product for complete Instructions, Warnings, Precautions, Possible Adverse Effects and Contraindications prior to use of the product.

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician.



AngioDynamics, the AngioDynamics logo, AlphaVac and the AlphaVac logo are trademarks and/or registered trademarks of AngioDynamics, Inc., an affiliate or subsidiary. ©2023 AngioDynamics, Inc. UC/VI/BR/1695 Rev 01 02/2023

